

CHEMSIL TR



PRIMJENA

tlak kemikalija - zbog nereagiranja na tjelesna tkiva i fluide, kao i mogućnosti sterilizacije parom, naročito se koriste u laboratorijima bolnica i ljekarni

MATERIJAL CIJEVI:

prozirni ekstrudirani silikon medicinske kvalitete

TEMPERATURA:

-40 +200 °C

STANDARD:

FDA 177.2600 (USA)

Promjer (mm)	Težina (kg/m)	Kolut (m)
3 x 5	0.015	100
3 x 6	0.025	100
4 x 6	0.019	100
4 x 7	0.031	100
4 x 8	0.045	100
5 x 8	0.037	100
6 x 9	0.042	100
8 x 12	0.058	100
10 x 15	0.072	100